



SURGI TRACTOR

■ *Innovative* ■ *Advanced* ■ *Wound Retractor / Protector*

 **SurgiCore**
Core in medical business

General

SurgiTractor designed to achieve the maximum opening with minimum incision maintaining proper moisture and protecting the contamination of the incised wound.


Features

Maximum opening with minimum incision

Thanks to the structure of upper ring of SurgiTractor, it gives the power to maintain the pure circle when it is placed on the operating site.

Easy to fold

Due to the Core Spring inside of the upper tube, the operator can easily fold the elastic film over the upper tube during the retracting stage.



Ordering Information

Ref. No.	Size	Description(Recommended incision size)	Contents	
			Inbox(ea)	Cartonbox(ea)
ST0103	XX-Small	SurgiTractor, disposable elastic wound retractor(1cm~3cm)	10	120
ST0204	X-Small	SurgiTractor, disposable elastic wound retractor(2cm~4cm)	10	120
ST0306	Small	SurgiTractor, disposable elastic wound retractor(3cm~6cm)	10	100
ST0509	Medium	SurgiTractor, disposable elastic wound retractor(5cm~9cm)	8	80
ST0914	Large	SurgiTractor, disposable elastic wound retractor(9cm~14cm)	8	48
ST1117	X-Large	SurgiTractor, disposable elastic wound retractor(11cm~17cm)	5	20
ST1725	XX-Large	SurgiTractor, disposable elastic wound retractor(17cm~25cm)	6	12
ST0306K	Small	SurgiTractor, Kit for a small incision(3cm~6cm), contents; 1 unit of small size wound retractor and a silicone cover	6	36
ST0306KA	Small	SurgiTractor, Kit for a small incision(3cm~6cm), contents; 1 unit of small size wound retractor and a silicone cover without a gas port	6	36
ST0306KB	Small	SurgiTractor, Kit for a small incision(3cm~6cm), contents; 1 unit of small size wound retractor and a silicone cover with a hook type cover-ring	6	36

Accessory

Ref. No.	Size	Description(Recommended incision size)	Contents	
			Inbox(ea)	Cartonbox(ea)
STC0306	Small	Silicone cover includes a sealed 12mm port and a Gas port, for use with ST0306	6	36

Manufacturer

SurgiCore Co., Ltd.
43-10, Jingoksandan 5beon-Ro, Gwangsan-Gu,
Gwangju 62207 Republic of Korea

Tel. +82-62-714-3850

Fax. +82-62-714-3861

www.surgicore.co.kr

Distributor

Versium, s.r.o
tel: +421 2 50 100 993,
info@versium.sk,
www.versium.sk

