

made in Germany

## Fast - Effective - Safe



180 mm Applicator



Flexible



Dimensionally stable

Atraumatic applicator tip

# HaemoCer<sup>™</sup> Flex SET



# BioCer

Flex SET Hoemostotic Powder + FlexApplicotor For all minimally invasive surgeries



## **Haemostatic Powder**

- Treatment of diffuse and localized bleeding
- For bleeding around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax in bone bleeding
- Reduction in the need for electrocautery and associated tissue damage
- Use with compatible filter cell-savers

## HaemoCer<sup>™</sup> FlexApplicator

- Longer applicator for bleeding sites that are difficult to access
- HaemoCer<sup>™</sup> FlexApplicator is bendable at several points at the same time and dimensionally stable to reach even difficult surgical sites
- The transparency of the applicator enables a clear view to the amount of HaemoCer<sup>™</sup> powder applied, therefore a controlled and pinpointed application is guaranteed
- Even after multiple bendings the powder freely flows through the applicator
- Atraumatic applicator tip prevents surgical lesions

### Product Overview HaemoCer<sup>™</sup> Flex Set

HFS215	HaemoCer™ 5g + 180 mm Applicator	1 pcs
HFS225	HaemoCer™ PLUS 5g + 180 mm Applicator	1 pcs

#### **BioCer Entwicklungs-GmbH** Ludwig-Thoma-Str. 36c • 9544

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Nemecko · Tel. +49 (0)921 78 77 70 0



#### **Distribútor pre SR:** Versium, s.r.o tel: +421 2 50 100 993, info@versium.sk, www.versium.sk

